

# Shot Noise

## Why is noise interesting?

In a conductivity measurement, noise can give a lot of **precious information** about:

- the **charge of (quasi-)particles** that contribute to the current (e.g. charge  $e/3$  in the fractional quantum Hall regime),
- the **statistics** obeyed by these particles (fermions, bosons),
- **interactions causing correlations** between them (e.g. Cooper pairs).


## The experiment

In this basic experiment, we propose you to determine the value of the elementary charge  $|e|$ .

## What is shot noise?



**Shot noise** arises because electronic **charge is transported in quantized portions**.

Example: electrons with a net charge  $-e$



## How can we measure it?

As in real lab work, you will be in charge of finding **good experimental conditions**.



Vacuum tube =  
noise generator =  
damped harmonic  
oscillator (RLC)

Voltage amplifier  
→ amplify  $U_{out}^{noise}$

Rectifier →  
Rectify  $U_{out}^{noise}$  into a  
DC voltage  $U_{amp}$

Voltmeter  
→ measure  $U_{amp}$

- Look for the  $f_{res}$  of this circuit!  
→ You will obtain  $C_{tot}$ ...
- Determine  $R$  with the help of the RLC bridge!

- You can read  $P...$
- And determine the amplification factor  $GA^2$  thanks to the HF oscillator and thus  $\langle U_{amp}^2(t) \rangle$ !

$$|e| = \frac{PC_{tot}}{GA^2R\langle I \rangle} \quad +/- \text{ some experimental error...}$$

